

INSTRUCTIONS FOR USE

Ref: APEX-IFU-CUP0300

Description	Specimen Cup 300ml Disposable Medical Pulp
Product Part Number	CUP0300
Ordering Item Code	101xxxCUP0300 (xxx – To be assigned upon order)
UDI-DI	95510204302438

Product Image



Intended Use

For patient treatment and monitoring, the Food Grade Specimen Cup with 300ml rated capacity is made from eco-friendly biodegradable sugarcane bagasse, intended for use as a multi-purpose cup. It is designed with smooth inner finishing and meshed outer surface for extra grip. Do **NOT** pour alcohol-based solutions into the Specimen Cup.


Storage

Keep in a cool, dry place, away from heat, flames, and sunlight.

User Guide



1
Collect urine sample or solution in the Specimen Cup.



2
After use, dispose the Specimen Cup into macerator or appropriate waste bins.



WARNING

Any serious incident related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user is established.

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.



APEX PULP Pte Ltd

Company Registration No.: 200715327H
132 Joo Seng Road #04-01. Singapore 368358
Tel: +65 6817 8969 | Fax: +65 6636 9813
www.apexpulp.com | sales@apexpulp.com

