

INSTRUCTIONS FOR USE

Ref: APEX-IFU-TIJ0350

Description	Injection Tray 350ml Disposable Medical Pulp
Product Part Number	TIJ0350
Ordering Item Code	101xxxTIJ0350 (xxx – To be assigned upon order)
UDI-DI	95510204302292

Product Image




Intended Use

For patient treatment and monitoring, the Food Grade Injection Tray with 350ml rated capacity is made from eco-friendly biodegradable sugarcane bagasse, intended for use as multi-purpose tray to carry syringes and other medical items. It is designed with smooth inner finishing and meshed outer surface for extra grip. Do **NOT** pour alcohol-based solutions into the Injection Tray.


Storage

Keep in a cool, dry place, away from heat, flames, and sunlight.

User Guide



Place surgical items on the Injection Tray.



After use, dispose the Injection Tray into macerator or appropriate waste bins.



WARNING

Any serious incident related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user is established.

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.



APEX PULP Pte Ltd

Company Registration No.: 200715327H
132 Joo Seng Road #04-01. Singapore 368358
Tel: +65 6817 8969 | Fax: +65 6636 9813
www.apexpulp.com | sales@apexpulp.com

