

INSTRUCTIONS FOR USE

Ref: APEX-IFU-GPB1000

Description	General Purpose Bowl 1000ml Disposable Medical Pulp
Product Part Number	GPB1000
Ordering Item Code	100xxxGPB1000 (xxx – To be assigned upon order)
UDI-DI	95510204300694

Product Image



Intended Use

For patient treatment and monitoring, the General Purpose Bowl Graduated with 1L rated capacity is made from eco-friendly biodegradable kraft material, intended for use as a general bowl for washing or vomiting. It is marked with graduation levels ranging from 200ml to 1000ml. Do **NOT** dispose towel into macerator.

Storage

Keep in a cool, dry place, away from heat, flames, and sunlight.


User Guide




Use General Purpose Bowl for washing or Vomiting.



After use, empty the contents into sink or toilet bowl.



Dispose the General Purpose Bowl directly into macerator.



Dispose the General Purpose Bowl directly into appropriate waste bins.



WARNING

Any serious incident related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user is established.

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.



APEX PULP Pte Ltd

Company Registration No.: 200715327H
132 Joo Seng Road #04-01. Singapore 368358
Tel: +65 6817 8969 | Fax: +65 6636 9813
www.apexpulp.com | sales@apexpulp.com

