

Date: 20/08/2022 Revision Date: 30/09/2025

Ref: APEX-IFU-GPB3000

Revision No.: 1

# INSTRUCTIONS FOR USE

Description **Product Part Number** Ordering Item Code **UDI-DI** 

General Purpose Bowl 3000ml Disposable Medical Pulp **GPB3000** 

100xxxGPB3000 (xxx – To be assigned upon order) 95510204300830

**Product Image** 















Intended Use

For patient treatment and monitoring, the General Purpose Bowl Graduated with 3L rated capacity is made from eco-friendly biodegradable kraft material, intended for use as a general bowl to rinse towels and/or to hold several litres of water. It is marked with graduation levels ranging from 500ml to 3000ml. Do NOT dispose towel into macerator.

#### Storage

Keep in a cool, dry place, away from heat, flames, and sunlight.

# User Guide




Fill the General Purpose Bowl with wate to the required level.



Wet or rinse the towel with water inside the bowl.



After use, pour away the water.



Dispose the General **Purpose Bowl directly** into macerator or into appropriate waste bins.



### WARNING

Any serious incident related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user is established.

## Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.











