

INSTRUCTIONS FOR USE

Ref: APEX-IFU-LMS2000

Description	Lid for Maxi Slipper Pan Liner Disposable Medical Pulp
Product Part Number	LMS2000
Ordering Item Code	100xxxLMS2000 (xxx – To be assigned upon order)
UDI-DI	<pending>

Product Image




Intended Use

Designed to cover the Maxi Slipper Pan Liner (MSP2000), helping prevent splashes and hide waste contents from sight and smell. Made from biodegradable paper material, it is macerator-friendly for convenient and eco-friendly disposal. Do **NOT** dispose of any plastic supports into the macerator.


Storage

Keep in a cool, dry place, away from heat, flames, and sunlight.

User Guide



After use, remove the Maxi Slipper Pan Liner from the plastic support, place a lid cover over it to contain the odour and sight.



Dispose the lid and Maxi Slipper Pan Liner directly into macerator or empty the contents from Maxi Slipper Pan Liner before disposing into an appropriate waste bin.



WARNING

Any serious incident related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user is established.

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.



APEX PULP Pte Ltd

Company Registration No.: 200715327H
132 Joo Seng Road #04-01. Singapore 368358
Tel: +65 6817 8969 | Fax: +65 6636 9813
www.apexpulp.com | sales@apexpulp.com

