

INSTRUCTIONS FOR USE

Ref: APEX-IFU-LMU1000

Description	Lid for Male Urinal Bulbous Disposable Medical Pulp
Product Part Number	LMU1000
Ordering Item Code	100xxxLMU1000 (xxx – To be assigned upon order)
UDI-DI	95510204303114


Product Image



Intended Use	For patient treatment and monitoring, the Lid for Male Urinal Bulbous is made from eco-friendly biodegradable kraft material, intended for use as a cap fitted over the opening of male urinal bulbous to prevent spillage and reduce odour.
--------------	--

Storage	Keep in a cool, dry place, away from heat, flames, and sunlight.
---------	--

User Guide



Remove the lid from Male Urinal Bulbous before use.



Place the Male Urinal Bulbous between the legs for urination.



Dispose the Male Urinal Bulbous with lid directly into macerator or empty the content from the Male Urinal Bulbous before disposing into appropriate waste bins.



WARNING

Any serious incident related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user is established.

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.



APEX PULP Pte Ltd

Company Registration No.: 200715327H
132 Joo Seng Road #04-01. Singapore 368358
Tel: +65 6817 8969 | Fax: +65 6636 9813
www.apexpulp.com | sales@apexpulp.com

