

INSTRUCTIONS FOR USE

Ref: APEX-IFU-MUB1000

Description	Male Urinal Bulbous 1000ml Disposable Medical Pulp
Product Part Number	MUB1000
Ordering Item Code	100xxxMUB1000 (xxx – To be assigned upon order)
UDI-DI	95510204301684

Product Image



Intended Use

For patient treatment and monitoring, the Male Urinal Bulbous with 1L rated capacity is made from eco-friendly biodegradable kraft material, intended for use as male urine bottle with large opening for easy access and long neck for patient to hold it comfortably. It is designed with uniform and sturdy base for it to stand on flat surfaces. It can also be fitted with an adaptor for use as Female Urinal Bulbous.

Storage

Keep in a cool, dry place, away from heat, flames, and sunlight.


User Guide



Place the Male Urinal Bulbous between the legs for urination.



Dispose the Male Urinal Bulbous directly into macerator or empty the content from the Male Urinal Bulbous before disposing into appropriate waste bins.



WARNING

Any serious incident related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user is established.

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.



APEX PULP Pte Ltd

Company Registration No.: 200715327H
132 Joo Seng Road #04-01. Singapore 368358
Tel: +65 6817 8969 | Fax: +65 6636 9813
www.apexpulp.com | sales@apexpulp.com

