

## GRADUATED MEASURING JUG 1000ML


**Eco-Friendly Medical Pulp** 



Designed for Single Use to Reduce Risk of Infection



## **Graduated Measuring Jug 1000ml Eco-Friendly Medical Pulp**



Manufactured as single-use disposable medical pulp to meet the requirement to hold water and/or human waste for a minimum period of four hours at an initial temperature of 35°C. Suitable for disposal into macerator which complies with PAS29:1999 "Disposable Pulp For Use in Healthcare".

## Suitable for







**Nursing Homes** 












## **User Guide**



Place the Graduated Measuring Jug on a flat surface, pour liquid into it slowly to prevent splashes.



Read and record volume of the liquid from the level markings.



Dispose the Graduated Measuring Jug directly into macerator or empty liquid from the jug before disposing into appropriate waste bins.

3b



Stackable



Water Retention



Macerator Friendly



**Eco-Friendly** 



Recycled

**PRODUCT DESCRIPTION** 

**PRODUCT SIZE** 

**PACKAGING** 

**CARTON DIMENSION** 

**ORDER CODE** 

Graduated Measuring Jug 1000ml

260 x 180 x 115 mm

100 pcs per carton

510 x 260 x 360 mm

100APPGMJ1000







Manufactured in accordance with PAS29:1999 Disposable Pulp Products For Use in Healthcare



