

LID FOR BEDPAN LINER & MIDI SLIPPER PAN LINER

Eco-Friendly Medical Pulp




Designed for Single Use to
Reduce Risk of Infection

Lid for Bedpan Liner and Midi Slipper Pan Liner Biodegradable Medical Pulp

This lid provides a hygienic way to cover our Bedpan Liners (BPG2000) and Midi Slipper Pan Liners (MSP1300), helping reduce odors after use. Fully macerator-friendly and made from biodegradable materials, it supports sustainable waste management practices in healthcare settings.

Manufactured in accordance with PAS29:1999 Disposable Pulp Products For Use in Healthcare.



Slits and cut-outs help to catch the macerator blades and makes it easier to tear and shred.




Lid with Bedpan Liner Graduated (BPG2000)




Lid with Midi Slipper Pan Liner (MSP1300)

User Guide



After use, remove the Bedpan Liner with support, place a lid cover over it to contain the odour and sight.



Dispose the Bedpan Liner with disposable support directly into macerator or empty the content from Bedpan Liner before disposing into appropriate waste bins.

Suitable for



Hospitals



Nursing Homes



Personal Use



Medical Clinics



Disability Homes

Note: DO NOT dispose the plastic support into macerator.



Single-Use



Macerator Friendly



Eco-Friendly



Water Retention



Stackable



Recycled



No Direct Sunlight



Keep Dry

PRODUCT DESCRIPTION

Lid for Bedpan Liner / Midi Slipper Pan Liner

PRODUCT SIZE

330 x 287 x 0.4mm

PACKAGING

200 pcs per carton

CARTON DIMENSION

340 x 300 x 90mm

UDI-DI: 95510204302988

ORDER CODE: 100xxxLBP2000 (xxx to be assigned upon ordering)



APEX PULP PTE LTD

Company Registration No.: 200715327H
132 Joo Seng Road #04-01, Singapore 368358
Tel: +65 6817 8969 | Fax: +65 6636 9813
www.apexpulp.com | sales@apexpulp.com

DOC. REF: APEX-PIB-LBP2000

Date of Revision: 28/01/2026

Revision: 0