

LID FOR MALE URINAL BULBOUS

Eco-Friendly Medical Pulp



Made from
BIODEGRADABLE
Double Lined Kraft

Designed for Single Use to Reduce Risk of Infection



Lid for Male Urinal Bulbous Eco-Friendly Medical Pulp





The Lid for Male Urinal Bulbous is made from eco-friendly biodegradable kraft material, intended for use as a cap fitted over the opening of male urinal bulbous to prevent spillage and reduce odour.

Manufactured as single-use disposable medical pulp to meet the requirement to hold water and/or human waste for a minimum period of four hours at an initial temperature of 35°C. Suitable for disposal into macerator which complies with PAS29:1999 "Disposable Pulp For Use in Healthcare".

Suitable for













User Guide



Remove the lid from Male Urinal Bulbous and before use.



Place the Male Urinal Bulbous between the legs for urination.





Dispose the Male Urinal Bulbous together with the lid directly into macerator or empty the content from the Male Urinal Bulbous before disposing into appropriate waste bins.









| PROD | UCT | DESC | RIPT | ION |
|------|-----|------|------|-----|
| | | | | |

PRODUCT SIZE

PACKAGING

CARTON DIMENSION

ORDER CODE

Lid for Male Urinal Bulbous

 $85 \times 71 \times 17 \text{ mm}$

1000 pcs per carton

380 x 340 x 260 mm

100APPLMU1000









Manufactured in accordance with PAS29:1999
Disposable Pulp Products For Use in Healthcare

