


ANAESTHETIC TRAY 2-COMPARTMENT

Eco-Friendly Medical Pulp



Designed for Single Use to Reduce Risk of Infection






User Guide



Place surgical items on the Anaesthetic Trau.



After use, dispose the Anaesthetic Tray into macerator or appropriate waste bins.

Anaesthetic Tray 2-Compartment Eco-Friendly Medical Pulp

The Food Grade Anaesthetic Tray with 2-Compartment is made from eco-friendly biodegradable sugarcane bagasse, intended for use as multi-purpose tray to carry ampoules at the small compartment, surgical instruments, dressings, wipes, creams and other medical items at the large compartment. It is designed with smooth inner finishing and meshed outer surface for extra grip.

Manufactured as single-use disposable medical pulp to meet the requirement to hold water and/or human waste for a minimum period of four hours at an initial temperature of 35°C. Suitable for disposal into macerator which complies with PAS29:1999 "Disposable Pulp For Use in Healthcare".

Features

- Consist of two deep compartments, ideal for surgical procedures
- Able to retain water for at least 4 hours
- Non-sterile, single use trau
- Meshed outer surface for better grip
- Made from sugarcane bagasse
- Biodegradable and compostable
- Can be disposed into macerator
- Manufactured in accordance with:
 - O EU Regulation No. 1935/2004 Materials intended to come into contact with food
 - O US FDA 21 CFR 176.170 Components of paper and paperboard in contact with aqueous and fatty foods.

Caution

Do not pour alcohol based solutions into the 2-Compartment Tray.











PRODUCI	DESCRIP	ION

PRODUCT SIZE

PACKAGING

CARTON DIMENSION

PRODUCT CODE

Anaesthetic Tray 2-Compartment (Food Grade)

225 x 180 x 50 mm

300 pcs per carton (6 x 50 pcs per bag) 510 x 260 x 210 mm

101APPT2C2000









Manufactured in accordance with PAS29:1999 Disposable Pulp Products For Use in Healthcare

