

## INSTRUCTIONS FOR USE

Ref: APEX-IFU-GWB3000

Description General Wash Bowl 3000ml Disposable Medical Pulp

Product Part Number GWB3000

Ordering Item Code 100APPGWB3000

Product Image



**Intended Use** For the purpose of patient treatment and monitoring, the General Purpose Bowl Graduated with 3L rated capacity is made from eco-friendly biodegradable kraft material, intended for use as a general bowl to rinse towels and/or to hold several litres of water. It is marked with graduation levels ranging from 500ml to 3000ml.

**Caution** Do not dispose the towel into macerator


### User Guide




1 Dispense and apply soap directly on the patient's body or on a towel.



2 Rinse the towel in the General Wash Bowl to wipe and cleanse patient's body.



3 After use, empty the contents into the sink.



4 Dispose the General Purpose Bowl directly into macerator or into appropriate waste bins.

### Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev2(1123)

