

INSTRUCTIONS FOR USE

Ref: APEX-IFU-LMU1000

Description Lid for Male Urinal Bulbous Disposable Medical Pulp

Product Part Number LMU1000

Ordering Item Code 100APPLMU1000


Product Image



Intended Use

For the purpose of patient treatment and monitoring, the Lid for Male Urinal Bulbous is made from eco-friendly biodegradable kraft material, intended for use as a cap fitted over the opening of male urinal bulbous to prevent spillage and reduce odour.

User Guide



1
Remove the lid from Male Urinal Bulbous before use.



2
Place the Male Urinal Bulbous between the legs for urination.



3a
Dispose the Male Urinal Bulbous with lid directly into macerator or empty the content from the Male Urinal Bulbous before disposing into appropriate waste bins.



Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev2(1123)

