

## INSTRUCTIONS FOR USE

Ref: APEX-IFU-MUB1000

Description Male Urinal Bulbous 1000ml Disposable Medical Pulp

Product Part Number MUB1000

Ordering Item Code 100APPMUB1000

Product Image



Intended Use

For the purpose of patient treatment and monitoring, the Male Urinal Bulbous with 1L rated capacity is made from eco-friendly biodegradable kraft material, intended for use as male urine bottle with large opening for easy access and long neck for patient to hold it comfortably. It is designed with uniform and sturdy base for it to stand on flat surfaces. It can also be fitted with an adaptor for use as Female Urinal Bulbous.


User Guide



Place the Male Urinal Bulbous between the legs for urination.



Dispose the Male Urinal Bulbous directly into macerator or empty the content from the Male Urinal Bulbous before disposing into appropriate waste bins.



### Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev2(1123)

