

INSTRUCTIONS FOR USE

Ref: APEX-IFU-TIJ0200

Description Injection Tray 200ml Disposable Medical Pulp

Product Part Number TIJ0200

Ordering Item Code 101APPTIJ0200

Product Image




Intended Use For the purpose of patient treatment and monitoring, the Food Grade Injection Tray with 200ml rated capacity is made from eco-friendly biodegradable sugarcane bagasse, intended for use as multi-purpose tray to carry syringes and other medical items. It is designed with smooth inner finishing and meshed outer surface for extra grip.

Caution Do not pour alcohol based solutions into the Injection Tray.

User Guide



1 Place surgical items on the Injection Tray.



2 After use, dispose the Injection Tray into macerator or appropriate waste bins.

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev2(1123)

