

INSTRUCTIONS FOR USE

Ref: APEX-IFU-CUP0300

Description Specimen Cup 300ml Disposable Medical Pulp

Item Code 101APPCUP0300

Product Image




Intended Use The Food Grade Specimen Cup with 300ml rated capacity is made from eco-friendly biodegradable sugarcane bagasse, intended for use as a multi-purpose cup. It is designed with smooth inner finishing and meshed outer surface for extra grip.

Caution Do not pour alcohol based solutions into the Specimen Cup.

User Guide



1 Collect urine sample or solution in the Specimen Cup.



2 After use, dispose the Specimen Cup into macerator or appropriate waste bins.

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev0(0822)

