

INSTRUCTIONS FOR USE

Ref: APEX-IFU-TIS1000

Description	Instrument Tray 1000ml Disposable Medical Pulp
Product Part Number	TIS1000
Ordering Item Code	101xxxTIS1000 (xxx – To be assigned upon order)
UDI-DI	95510204302124


Product Image




Intended Use	For patient treatment and monitoring, the Food Grade Instrument Tray with 1000ml rated capacity is made from eco-friendly biodegradable sugarcane bagasse, intended for use as multi-purpose tray to carry surgical instruments, dressings, wipes, creams and other medical items. It is designed with smooth inner finishing and meshed outer surface for extra grip. Do NOT pour alcohol-based solutions into the tray.
--------------	--

Storage	Keep in a cool, dry place, away from heat, flames, and sunlight.
---------	--

User Guide



Place surgical items on the Instrument Tray.



After use, dispose the Instrument Tray into macerator or appropriate waste bins.



WARNING

Any serious incident related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user is established.

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.



APEX PULP Pte Ltd

Company Registration No.: 200715327H
132 Joo Seng Road #04-01. Singapore 368358
Tel: +65 6817 8969 | Fax: +65 6636 9813
www.apexpulp.com | sales@apexpulp.com

