

BEDPAN SUPPORT <100 KG

Eco-Friendly Medical Pulp



Made from
BIODEGRADABLE
Double Lined Kraft

**Designed for Single Use to
Reduce Risk of Infection**

Bedpan Support <100kg Biodegradable Medical Pulp

The sturdy Bedpan Support that supports up to 100kg evenly distributed weight, is made from eco-friendly biodegradable kraft material, intended for use as support with bedpan liner place over it. Due to the innovative mould pulp design structure, it can withstand patient weighing up to 100kg when placed under the buttock.

Manufactured in accordance with PAS29:1999 Disposable Pulp Products For Use in Healthcare.

Suitable for



Hospitals



Nursing Homes



Personal Use




Medical Clinics




Disability Homes


User Guide




For user weight under 50kg, place bedpan liner over the disposable moulded pulp bedpan support.



If the user weight is between 50 to 100kg, stack 2 pieces together for added strength.



Slide the bedpan liner with support under the buttock. To prevent spills, place a sheet of incontinence pad under the bedpan.



Dispose the Bedpan Liner with disposable support directly into macerator or empty the content from Bedpan Liner before disposing into appropriate waste bins.



Single-Use



Macerator Friendly



Eco-Friendly



Water Retention



Stackable



Recycled



No Direct Sunlight



Keep Dry

PRODUCT DESCRIPTION	PRODUCT SIZE	PACKAGING	CARTON DIMENSION
Bedpan Support < 100kg	370 x 295 x 95mm	40 pcs per carton	755 x 370 x 310mm

UDI-DI: 95510204300144
ORDER CODE: 100xxxBPS2000 (xxx to be assigned upon ordering)

