



# APEX PULP

Manufacturer of Moulded Pulp



Made from  
**BIODEGRADABLE**  
Double Lined Kraft

<p><b>Male Urinal Bulbous 1000ml</b></p> <p><b>MUB1000</b></p> <p>100</p> 	<p><b>Adaptor for Urinal Bulbous (Female Urinal)</b></p> <p><b>AUB1000</b></p> <p>200</p> 	<p><b>Lid for Male Urinal Bulbous</b></p> <p><b>LMU1000</b></p> <p>1000</p> 	<p><b>Graduated Measuring Jug 1000ml</b></p> <p><b>GMJ1000</b></p> <p>100</p> 
<p><b>Kidney Dish 700ml</b></p> <p><b>KND0700</b></p> <p>300</p> 	<p><b>Kidney Dish 950ml</b></p> <p><b>KND0950</b></p> <p>300</p> 	<p><b>General Purpose Bowl 1000ml</b></p> <p><b>GPB1000</b></p> <p>100</p> 	<p><b>General Semi Bowl 1500ml</b></p> <p><b>GSB1500</b></p> <p>100</p> 
<p><b>General Purpose Pot 2000ml</b></p> <p><b>GPP2000</b></p> <p>100</p> 	<p><b>General Purpose Bowl 3000ml</b></p> <p><b>GPB3000</b></p> <p>100</p> 	<p><b>General Wash Bowl 3000ml</b></p> <p><b>GWB3000</b></p> <p>100</p> 	<p><b>General Wash Bowl 4000ml</b></p> <p><b>GWB4000</b></p> <p>100</p> 
<p><b>Bedpan Liner Graduated 2000ml</b></p> <p><b>BPG2000</b></p> <p>100</p> 	<p><b>Bedpan Support &lt;100kg</b></p> <p><b>BPS2000</b></p> <p>45</p> 	<p><b>Maxi Slipper Pan Liner 2000ml</b></p> <p><b>MSP2000</b></p> <p>100</p> 	<p><b>Midi Slipper Pan Liner 1300ml</b></p> <p><b>MSP1300</b></p> <p>100</p> 





Made from  
**BIODEGRADABLE**  
Sugarcane Bagasse

<b>Kidney Dish 950ml</b>	<b>Anaesthetic Tray 2-Compartment</b>	<b>Anaesthetic Tray 4-Compartment</b>	<b>Instrument Tray 1000ml</b>
<b>KNS0950</b> 300	<b>T2C2000</b> 300	<b>T4C2000</b> 100	<b>TIS1000</b> 500
<b>Injection Tray 200ml</b>	<b>Injection Tray 350ml</b>	<b>Medicine Cup 50ml</b>	<b>Gallipot 90ml</b>
<b>TIJ0200</b> 1000	<b>TIJ0350</b> 1000	<b>CUP0050</b> 1000	<b>CUP0090</b> 1000

### Features

- Able to retain water for at least 4 hours
- Non-sterile, single use
- Biodegradable and compostable
- Disposable into macerator or appropriate waste bin
- Manufactured in accordance with PAS29:1999  
Disposable Pulp Products For Use in Healthcare

Description & Features →

Item Code →

Quantity per carton →

	<b>Specimen Cup 300ml</b>
<b>CUP0300</b>	
<b>1000</b>	

### Suitable for





APEX PULP PTE. LTD.  
132 Joo Seng Road #04-01. Singapore 368358.  
Tel: +65 6817 8969      Fax: +65 6636 9813  
www.apexpulp.com      sales@apexpulp.com

**Manufactured in Malaysia**



+65 9820 4940