

INSTRUCTIONS FOR USE

Ref: APEX-IFU-TIS1000

Description Instrument Tray 1000ml Disposable Medical Pulp

Product Part Number TIS1000

Ordering Item Code 101APPTIS1000

Product Image




Intended Use

For the purpose of patient treatment and monitoring, the Food Grade Instrument Tray with 1000ml rated capacity is made from eco-friendly biodegradable sugarcane bagasse, intended for use as multi-purpose tray to carry surgical instruments, dressings, wipes, creams and other medical items. It is designed with smooth inner finishing and meshed outer surface for extra grip.


Caution

Do not pour alcohol based solutions into the Instrument Tray.

User Guide



Place surgical items on the Instrument Tray.



After use, dispose the Instrument Tray into macerator or appropriate waste bins.

Disclaimer:

All products specifications, illustrations, drawings, images, dimensions, performance data and other information on our Instructions for Use or made available by us are subjected to change without notice and they are intended to represent no more than a general illustration of the products intended use and their features and do not constitute a warranty or representation by us that the products will conform with the same. Products that are shipped to you may vary from these images or specifications.

Rev2(1123)

